

Orthoplex White Label MultiGen Biotic Stability

MultiGen Biotic capsules and powder have been formulated for clinical practitioners only and has been designed to deliver exceptional quality by ensuring a guaranteed 30 Billion Colony forming Units (CFU's) in each capsules or serve of powder. To preserve this potency and efficacy by the time of product expiry it is recommended **MultiGen Biotic** be kept refrigerated below 8°C.

However, MultiGen Biotic has been specifically formulated to be tolerant to short term temperature increases to allow for transit outside of refrigeration. Stability testing has shown that the CFU count in MultiGen Biotic is above the stated levels on the label after 9 days at a temperature up to 30°C and 2 days at 40°C.

So, whilst transporting with icepacks is desirable, MultiGen Biotic **sent in unrefrigerated packs, or left on the shelf, has the same quantity of organisms for up to 9 days at 30°C, and displays no discrepancies with efficacy or potency of the product.**

